

Product Datasheet

Histidine (6xHis-tag) Antibody HM1160

Product Name	Histidine (6xHis-tag) Antibody
Catalog Number	HM1160
Description	Purified Polyclonal Antibody
Species/Host	Mouse
Tested Applications	ELISA, LF
Immunogen	6xHis synthetic peptide
Brand	EastCoast Bio
Form/Appearance	Ion Exchange Purified Monoclonal Antibody
Storage	Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Buffer	0.015M Potassium Phosphate Buffer, pH 7.2 with 0.85% Sodium Chloride
Preservatives	0.05% Sodium Azide