

Product Datasheet

Influenza A Virus Antibody HM1163

Product Name	Influenza A Virus Antibody
Catalog Number	HM1163
Description	Purified Monoclonal Antibody
Species/Host	CHO Cells
Tested Applications	LF
Immunogen	Inactivated Influenza A Virus
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	+2-8 °C
Isotype	IgG2a
Clonality	Monoclonal
Purity	> 90 %
Buffer	10 mM Phosphate Buffered Saline, pH 7.4
Cross Reactivity	No reactivity with Influenza B, Adenovirus, RSV, Parainfluenza, M. pneumoniae or SARS-CoV-2
Specificity	Influenza A Virus Nucleoprotein
Matched Pairs	HM1162 (capture or detection)
Preservatives	0.1% Proclin 300