

## Product Datasheet

### Influenza B Virus Antibody HM1168

<b>Product Name</b>	Influenza B Virus Antibody
<b>Catalog Number</b>	HM1168
<b>Description</b>	Purified Recombinant Monoclonal Antibody
<b>Species/Host</b>	CHO Cells
<b>Tested Applications</b>	LF
<b>Immunogen</b>	Inactivated Influenza B Virus
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Purified Monoclonal Antibody
<b>Storage</b>	+2-8 °C
<b>Isotype</b>	IgG2b
<b>Clonality</b>	Monoclonal
<b>Purity</b>	> 90 %
<b>Buffer</b>	10 mM Phosphate Buffered Saline, pH 7.4
<b>Cross Reactivity</b>	No reactivity with Influenza A, Adenovirus, RSV, Parainfluenza, M. pneumoniae or SARS-CoV-2
<b>Specificity</b>	Influenza B Virus Nucleoprotein
<b>Matched Pairs</b>	HM1167 (capture)
<b>Preservatives</b>	0.1% Proclin 300