

# Product Datasheet

## Human CD106 ELISA Kit 850.580.096

Product Name	Human CD106 ELISA Kit
Catalog Number	850.580.096
Description	Diaclone ELISA Kits include pre-coated strip plates, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents, and a detailed procedure.
Range	1.56-50 ng/ml
Sensitivity	0.6 ng/ml
Tested Applications	ELISA
Alternative Names	VCAM-1
Brand	Diaclone
Storage	Store at +4°C
Cross Reactivity	No cross reactivity with other human adhesion molecules
Specificity	Recognizes both natural and recombinant human soluble VCAM-1
Sample Types	Serum, Plasma, Cell culture supernatant