

Product Datasheet

Mycobacteria (LAM) Antibody HM1177

Product Name	Mycobacteria (LAM) Antibody
Catalog Number	HM1177
Species/Host	Rabbit
Tested Applications	ELISA
Brand	EastCoast Bio
Form/Appearance	Purified Recombinant Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG
Clonality	Monoclonal
Purity	>90%
Buffer	10 mM Phosphate Buffered Saline, pH 7.2
Cross Reactivity	Reactivity with M. tuberculosis, M. smegmatis, M. avium, M. bovis, M. leprae, M. phlei. No reactivity with C. pneumoniae, M. pneumoniae, S.pneumoniae, L. pneumophila, K. pneumoniae, S. pyogenes
Specificity	Mycobacteria lipoarabinomannan (LAM)
Matched Pairs	HM1176 (capture)
Preservatives	0.1% Sodium Azide