

Product Datasheet

Adenovirus Antibody HM1178

Product Name	Adenovirus Antibody
Catalog Number	HM1178
Description	Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	WB, ELISA, LF, IFA
Brand	EastCoast Bio
Form/Appearance	Protein G Purified Monoclonal Antibody
Storage	Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Isotype	IgG2a
Clonality	Monoclonal
Buffer	25 mM Phosphate Buffered Saline, pH 7.0-7.5
Specificity	Hexon Protein
Matched Pairs	FG300 (capture)
Preservatives	0.1% Sodium Azide