

Product Datasheet

Human Chorionic Gonadotropin β (HCG β) Antibody HM118

Product Name	Human Chorionic Gonadotropin β (HCG β) Antibody
Catalog Number	HM118
Description	Protein A Purified Recombinant Monoclonal Antibody
Species/Host	Mouse
Tested Applications	ELISA
Immunogen	Human Chorionic Gonadotropin (beta subunit)
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	+2-8 °C
Isotype	IgG
Clonality	Monoclonal
Purity	> 95 %
Buffer	0.05M Phosphate Buffered Saline, pH 7.2 \pm 0.2
Cross Reactivity	β hCG 100%, ihCG 115%, LH < 1%, FSH < 1%, TSH < 1%
Specificity	HCG β subunit
Matched Pairs	HM119 (detection)
Preservatives	0.1% Sodium Azide