

Product Datasheet

Acrosomal Protein SP-10 (ACRV1) Antibody HM1188

Product Name	Acrosomal Protein SP-10 (ACRV1) Antibody
Catalog Number	HM1188
Description	Ion Exchange Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	LF, ELISA
Immunogen	Native Human Acrosomal Protein SP-10
Brand	EastCoast Bio
Form/Appearance	Ion Exchange Purified Monoclonal Antibody
Storage	Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Isotype	IgG2a
Clonality	Monoclonal
Purity	> 90 %
Buffer	0.015M Potassium Phosphate, pH 7.2 with 0.85% NaCl
Specificity	Human Acrosomal Protein SP-10
Preservatives	0.05% Sodium Azide