

## Product Datasheet

### Human Chorionic Gonadotropin $\beta$ (HCG $\beta$ ) Antibody HM119

<b>Product Name</b>	Human Chorionic Gonadotropin $\beta$ (HCG $\beta$ ) Antibody
<b>Catalog Number</b>	HM119
<b>Description</b>	Ion Exchange Purified Monoclonal Antibody
<b>Species/Host</b>	Mouse
<b>Tested Applications</b>	ELISA
<b>Immunogen</b>	Human Chorionic Gonadotropin (beta subunit)
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Purified Monoclonal Antibody
<b>Storage</b>	+2-8 °C
<b>Isotype</b>	IgG
<b>Clonality</b>	Monoclonal
<b>Purity</b>	> 95 %
<b>Buffer</b>	0.05 M Phosphate Buffered Saline, pH 7.4
<b>Cross Reactivity</b>	$\beta$ hCG 100%, ihCG 115%, LH < 1%, FSH < 1%, TSH < 1%
<b>Specificity</b>	HCG $\beta$ subunit
<b>Matched Pairs</b>	HM118 (capture)
<b>Preservatives</b>	0.1% Sodium Azide