

Product Datasheet

Human CD86 ELISA Kit 850.590.048

Product Name	Human CD86 ELISA Kit
Catalog Number	850.590.048
Description	Diaclone ELISA Kits include pre-coated strip plates, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents, and a detailed procedure.
Range	1.187-38 U/ml
Sensitivity	0.6 U/ml
Tested Applications	ELISA
Alternative Names	B7-2
Brand	Diaclone
Storage	Store at +4°C
Cross Reactivity	No cross reactivity with other human soluble molecules
Specificity	Recognizes natural human soluble CD86
Sample Types	Serum, Plasma, Cell culture supernatant