

Product Datasheet

NT-proBNP (38-44) Antibody HM1196

Product Name	NT-proBNP (38-44) Antibody
Catalog Number	HM1196
Description	Purified Monoclonal Antibody
Tested Applications	CLIA, LF
Immunogen	Native Human ProBNP antigen
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	-20 °C ± 5 °C
Isotype	IgG1
Clonality	Monoclonal
Purity	> 90 %
Buffer	10 mM Phosphate Buffer, pH 7.4 ± 0.2 with 0.15M NaCl
Specificity	NT-proBNP
Matched Pairs	HM148 (detection)
Preservatives	0.1% Proclin 300