

Product Datasheet

Human Immunodeficiency Virus 1 (HIV 1-p24) Antibody HM1202

Product Name	Human Immunodeficiency Virus 1 (HIV 1-p24) Antibody
Catalog Number	HM1202
Description	Protein A Purified Monoclonal Antibody
Tested Applications	ELISA, CLIA
Brand	EastCoast Bio
Form/Appearance	Monoclonal Antibody
Storage	-20 °C. Avoid repeated freezing and thawing
Clonality	Monoclonal
Purity	≥ 90% by SDS-PAGE
Buffer	10 mM phosphate buffer, pH 7.4 ± 0.2, with 0.15 M NaCl, 0.1% ProClin 300 as a preservative.
Specificity	HIV-1 p24