

# Product Datasheet

## C-Reactive Protein (CRP) Antibody HM121

Product Name	C-Reactive Protein (CRP) Antibody
Catalog Number	HM121
Species/Host	Mouse
Tested Applications	LF, ELISA
Immunogen	Purified human CRP
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term (≤ 2 weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Purity	>95%
Buffer	10 mM Phosphate Buffered Saline, pH 7.0
Specificity	Human CRP
Matched Pairs	HM122
Preservatives	0.09% Sodium Azide