

# Product Datasheet

## C-Reactive Protein (CRP) Antibody HM122

Product Name	C-Reactive Protein (CRP) Antibody
Catalog Number	HM122
Description	Protein A Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	LF, ELISA
Immunogen	Purified human CRP
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term (≤ 2 weeks): +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Purity	> 95 %
Buffer	10 mM Phosphate Buffered Saline, pH 7.0
Specificity	Human CRP
Matched Pairs	HM121
Preservatives	0.09% Sodium Azide