

Product Datasheet

Vitamin D Antibody HM127

Product Name	Vitamin D Antibody
Catalog Number	HM127
Species/Host	Mouse
Tested Applications	ELISA
Immunogen	Vitamin D-BSA
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term (≤ 2 weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Clonality	Monoclonal
Purity	>90%
Buffer	10 mM Phosphate Buffer pH 7.4 with 150 mM NaCl
Cross Reactivity	25OH-Vitamin D2 84%, 1,25(OH) ₂ -Vitamin D3 21%, 1,25(OH) ₂ -Vitamin D2 2.2%, Vitamin D3 2.9%, Vitamin D2 1.5%
Specificity	25OH-Vitamin D3
Preservatives	0.1% Proclin