

Product Datasheet

D-dimer Antibody HM135

Product Name	D-dimer Antibody
Catalog Number	HM135
Species/Host	Mouse
Tested Applications	ELISA, WB
Brand	EastCoast Bio
Form/Appearance	Protein A Purified Monoclonal Antibody
Storage	2-8°C
Isotype	IgG1
Clonality	Monoclonal
Buffer	10 mM Phosphate Buffered Saline, pH 7.4
Cross Reactivity	Does not cross-react with fibrinogen
Specificity	D-dimer, D-dimer β -chain, and fibrin degradation products
Matched Pairs	HM136 (detection)
Preservatives	0.09% Sodium Azide