

Product Datasheet

Human p53 ELISA Kit 850.630

Product Name	Human p53 ELISA Kit
Catalog Number	850.630
Description	Diaclone ELISA Kits include pre-coated strip plates, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents, and a detailed procedure.
Range	3.1-100 U/ml
Sensitivity	1.5 U/ml
Tested Applications	ELISA
Brand	Diaclone
Storage	Store at +4°C
Cross Reactivity	No cross reactivity with other human soluble molecules
Specificity	Recognizes both natural and recombinant soluble human p53
Sample Types	Serum, Plasma, Cell lysate