

Product Datasheet

Morphine Antibody HM139

Product Name	Morphine Antibody
Catalog Number	HM139
Description	Protein A Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	ELISA, LF, IFA
Brand	EastCoast Bio
Form/Appearance	Protein A Purified Monoclonal Antibody
Storage	Short Term (≤ 2 weeks): +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Isotype	IgG1, lamda
Clonality	Monoclonal
Purity	> 98 % by SDS-PAGE
Buffer	10 mM Phosphate Buffered Saline, pH 7.2 \pm 0.2
Specificity	Morphine
Preservatives	0.05% Sodium Azide