

Product Datasheet

Carbamazepine Antibody HM140

Product Name	Carbamazepine Antibody
Catalog Number	HM140
Description	Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	LF, ELISA
Brand	EastCoast Bio
Form/Appearance	Protein A Purified Monoclonal Antibody
Storage	Short Term (≤ 2 weeks): +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Isotype	IgG2a, kappa
Clonality	Monoclonal
Purity	> 98 % by SDS-PAGE
Buffer	10 mM Phosphate Buffered Saline, pH 7.2 \pm 0.2
Preservatives	0.05% Sodium Azide