

## Product Datasheet

### Gentamicin Antibody HM141

<b>Product Name</b>	Gentamicin Antibody
<b>Catalog Number</b>	HM141
<b>Description</b>	Purified Monoclonal Antibody
<b>Species/Host</b>	Mouse
<b>Tested Applications</b>	LF, ELISA, IFA
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Protein A Purified Monoclonal Antibody
<b>Storage</b>	Short Term ( $\leq 2$ weeks): +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
<b>Isotype</b>	IgG1, kappa
<b>Clonality</b>	Monoclonal
<b>Purity</b>	$\geq 98\%$ by SDS-PAGE
<b>Buffer</b>	Phosphate Buffered Saline, pH 7.2 $\pm$ 0.2
<b>Preservatives</b>	0.05% Sodium Azide