

Product Datasheet

Gentamicin Antibody HM141

Product Name	Gentamicin Antibody
Catalog Number	HM141
Description	Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	ELISA, LF, IFA
Brand	EastCoast Bio
Form/Appearance	Protein A Purified Monoclonal Antibody
Storage	Short Term (≤ 2 weeks): +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Isotype	IgG1, kappa
Clonality	Monoclonal
Purity	$\geq 98\%$ by SDS-PAGE
Buffer	Phosphate Buffered Saline, pH 7.2 \pm 0.2
Preservatives	0.05% Sodium Azide