

Product Datasheet

Creatine Kinase MB (CK-MB) Antibody HM144

Product Name	Creatine Kinase MB (CK-MB) Antibody
Catalog Number	HM144
Description	Protein A Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	LF, ELISA, IFA
Brand	EastCoast Bio
Form/Appearance	Protein A Purified Monoclonal Antibody
Storage	Short Term (≤ 2 weeks): +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Isotype	IgG2b, kappa
Clonality	Monoclonal
Purity	$\geq 98\%$ by SDS-PAGE
Buffer	10 mM Phosphate Buffered Saline, pH 7.2 \pm 0.2
Specificity	Creatine Kinase MB
Matched Pairs	HM1128 (detection)
Preservatives	0.05% Sodium Azide