

## Product Datasheet

### Creatine Kinase MB (CK-MB) Antibody HM144

<b>Product Name</b>	Creatine Kinase MB (CK-MB) Antibody
<b>Catalog Number</b>	HM144
<b>Description</b>	Protein A Purified Monoclonal Antibody
<b>Species/Host</b>	Mouse
<b>Tested Applications</b>	ELISA, LF, IFA
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Protein A Purified Monoclonal Antibody
<b>Storage</b>	Short Term ( $\leq 2$ weeks): +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
<b>Isotype</b>	IgG2b, kappa
<b>Clonality</b>	Monoclonal
<b>Purity</b>	$\geq 98\%$ by SDS-PAGE
<b>Buffer</b>	10 mM Phosphate Buffered Saline, pH 7.2 $\pm$ 0.2
<b>Specificity</b>	Creatine Kinase MB
<b>Matched Pairs</b>	HM1128 (detection)
<b>Preservatives</b>	0.05% Sodium Azide