

Product Datasheet

NT-proBNP (27-31) Antibody HM148

| | |
|----------------------------|--|
| Product Name | NT-proBNP (27-31) Antibody |
| Catalog Number | HM148 |
| Description | Purified Monoclonal Antibody |
| Species/Host | Mouse |
| Tested Applications | LF, CLIA |
| Immunogen | Native Human ProBNP antigen |
| Brand | EastCoast Bio |
| Form/Appearance | Purified Monoclonal Antibody |
| Storage | Short Term (≤ 2 weeks): +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing. |
| Isotype | IgG1 |
| Clonality | Monoclonal |
| Purity | > 90 % |
| Buffer | 10 mM Phosphate Buffer, pH 7.4 \pm 0.2 with 0.15M NaCl |
| Specificity | NT-proBNP |
| Matched Pairs | HM146, HM899 or HM1196 (capture) |