

Product Datasheet

NT-proBNP (27-31) Antibody HM148

Product Name	NT-proBNP (27-31) Antibody
Catalog Number	HM148
Description	Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	CLIA, LF
Immunogen	Native Human ProBNP antigen
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term (≤ 2 weeks): +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Purity	> 90 %
Buffer	10 mM Phosphate Buffer, pH 7.4 \pm 0.2 with 0.15M NaCl
Specificity	NT-proBNP
Matched Pairs	HM146, HM899 or HM1196 (capture)