

# Product Datasheet

## NT-proBNP (27-31) Antibody HM148

Product Name	NT-proBNP (27-31) Antibody
Catalog Number	HM148
Species/Host	Mouse
Tested Applications	LF, CLIA
Immunogen	Native Human ProBNP antigen
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term ( $\leq$ 2 weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Purity	>90%
Buffer	10mM Phosphate Buffered Saline, pH 7.4
Specificity	NT-proBNP
Matched Pairs	HM146, HM899 or HM1196 (capture)