

Product Datasheet

Human CD62L ELISA Kit 850.650.048

Product Name	Human CD62L ELISA Kit
Catalog Number	850.650.048
Description	Diaclone ELISA Kits include pre-coated strip plates, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents, and a detailed procedure.
Range	1.9-60 ng/ml
Sensitivity	1 ng/ml
Tested Applications	ELISA
Alternative Names	L-selectin / LAM-1
Brand	Diaclone
Storage	Store at +4°C
Cross Reactivity	No cross reactivity with other human adhesion molecules
Specificity	Recognizes both natural and recombinant human soluble LAM-1
Sample Types	Serum, Plasma, Cell culture supernatant