

Product Datasheet

Heartworm, Canine Antibody HM163

Product Name	Heartworm, Canine Antibody
Catalog Number	HM163
Species/Host	Mouse
Tested Applications	ELISA
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Purity	>90%
Buffer	10 mM Phosphate Buffered Saline, pH 7.2
Cross Reactivity	No reactivity with uninfected canine serum
Specificity	Dirofilaria immitis
Matched Pairs	HM162 or HM761 (detection)
Preservatives	0.1% Sodium Azide