

Product Datasheet

Human CD62L ELISA Kit 850.650.192

| | |
|----------------------------|--|
| Product Name | Human CD62L ELISA Kit |
| Catalog Number | 850.650.192 |
| Description | Diaclone ELISA Kits include pre-coated strip plates, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents, and a detailed procedure. |
| Range | 1.9-60 ng/ml |
| Sensitivity | 1 ng/ml |
| Tested Applications | ELISA |
| Alternative Names | L-selectin / LAM-1 |
| Brand | Diaclone |
| Storage | Store at +4°C |
| Cross Reactivity | No cross reactivity with other human adhesion molecules |
| Specificity | Recognizes both natural and recombinant human soluble LAM-1 |
| Sample Types | Serum, Plasma, Cell culture supernatant |