

# Product Datasheet

## Human CD116 ELISA Kit 850.660.192

Product Name	Human CD116 ELISA Kit
Catalog Number	850.660.192
Description	Diaclone ELISA Kits include pre-coated strip plates, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents, and a detailed procedure.
Range	31.25-1000 pg/ml
Sensitivity	16 pg/ml
Tested Applications	ELISA
Alternative Names	GM-CSF R
Brand	Diaclone
Storage	Store at +4°C
Cross Reactivity	No cross reactivity with other human soluble receptors
Specificity	Recognizes both natural and recombinant human soluble GM-CSF receptor a chain
Sample Types	Serum, Plasma, Cell culture supernatant