

Product Datasheet

Feline Leukemia Virus (p27) Antibody HM189

Product Name	Feline Leukemia Virus (p27) Antibody
Catalog Number	HM189
Description	Protein A Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	WB, ELISA, LF
Brand	EastCoast Bio
Form/Appearance	Protein A Purified Monoclonal Antibody
Storage	Short Term (≤ 2 weeks): +2-8 °C. Long Term: -20 °C or below. Avoid repeated freezing and thawing.
Isotype	IgG1, kappa
Clonality	Monoclonal
Buffer	10 mM Phosphate Buffer, pH 7.2
Specificity	FeLV p27
Matched Pairs	HM190 (capture)
Preservatives	0.1% Sodium Azide