

Product Datasheet

Troponin I (41-49) Antibody HM192

Product Name	Troponin I (41-49) Antibody
Catalog Number	HM192
Species/Host	Mouse
Tested Applications	LF, ELISA, WB
Brand	EastCoast Bio
Form/Appearance	Protein A Purified Monoclonal Antibody
Storage	2-8°C
Isotype	IgG2b
Clonality	Monoclonal
Buffer	10 mM Phosphate Buffered Saline, pH 7.4
Cross Reactivity	Does not cross-react with skeletal muscle Troponin I
Specificity	Troponin I (epitope 41-49 a.a.)
Preservatives	0.09% Sodium Azide