

Product Datasheet

Human CD124 ELISA Kit 850.670.048

Product Name	Human CD124 ELISA Kit
Catalog Number	850.670.048
Description	Diaclone ELISA Kits include pre-coated strip plates, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents, and a detailed procedure.
Range	31.25-1000 pg/ml
Sensitivity	11.5 pg/ml
Tested Applications	ELISA
Alternative Names	IL-4R
Brand	Diaclone
Storage	Store at +4°C
Cross Reactivity	No cross reactivity with other human soluble receptors
Specificity	Recognizes both natural and recombinant human soluble IL-4R
Sample Types	Serum, Plasma, Cell culture supernatant