

Product Datasheet

Distemper Virus, Canine Antibody HM196

Product Name	Distemper Virus, Canine Antibody
Catalog Number	HM196
Description	Protein A Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	ELISA, LF
Brand	EastCoast Bio
Form/Appearance	Protein A Purified Monoclonal Antibody
Storage	+2-8 °C
Isotype	IgG2a
Clonality	Monoclonal
Buffer	10 mM Phosphate Buffered Saline
Specificity	CDV surface envelope antigen
Matched Pairs	HM195 (detection)
Preservatives	0.09% Sodium Azide