

Product Datasheet

Listeria monocytogenes Antibody HM217

Product Name	Listeria monocytogenes Antibody
Catalog Number	HM217
Species/Host	Mouse
Tested Applications	ELISA, IFA
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG2b
Clonality	Monoclonal
Purity	>90%
Buffer	10 mM Phosphate Buffered Saline, pH 7.2
Specificity	L. monocytogenes p60
Matched Pairs	HM218 or HM219 (detection)
Preservatives	0.1% Sodium Azide