

Product Datasheet

Procalcitonin (PCT-katacalcin) Antibody HM220

Product Name	Procalcitonin (PCT-katacalcin) Antibody
Catalog Number	HM220
Description	Protein A Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	LF
Immunogen	Purified recombinant human Procalcitonin
Brand	EastCoast Bio
Form/Appearance	Protein A Purified Monoclonal Antibody
Storage	Short Term (≤ 2 weeks): +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Purity	> 95 % by SDS-PAGE
Buffer	10 mM Phosphate Buffered Saline, pH 7.2
Specificity	Katacalcin (C-terminal) epitope
Matched Pairs	HM078 (capture)
Preservatives	0.09% Sodium Azide