

Product Datasheet

Luteinizing Hormone α (LH- α) Antibody HM222

Product Name	Luteinizing Hormone α (LH- α) Antibody
Catalog Number	HM222
Description	Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	ELISA
Alternative Names	LH, Lutropin, Interstitial Cell-Stimulating Hormone
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	+2-8 °C
Isotype	IgG
Clonality	Monoclonal
Purity	> 95 %
Buffer	0.05M Phosphate Buffered Saline, pH 7.2 \pm 0.2
Cross Reactivity	LH 100%, β hCG < 1%, FSH < 1%, ihCG 7%, TSH < 1%
Specificity	LH α
Matched Pairs	FA313 or FA309 (capture)
Preservatives	0.1% Sodium Azide