

Product Datasheet

Human Serum Albumin Antibody HM223

Product Name	Human Serum Albumin Antibody
Catalog Number	HM223
Species/Host	Mouse
Tested Applications	ELISA, WB
Brand	EastCoast Bio
Form/Appearance	Protein A Purified Monoclonal Antibody
Storage	2-8°C
Isotype	IgG1
Clonality	Monoclonal
Buffer	10 mM Phosphate Buffered Saline, pH 7.4
Cross Reactivity	No reactivity with BSA, egg albumin, or human alpha fetoprotein
Specificity	Human Serum Albumin
Preservatives	0.09% Sodium Azide