

Product Datasheet

Human CD117 ELISA Kit 850.700.048

Product Name	Human CD117 ELISA Kit
Catalog Number	850.700.048
Description	Diaclone ELISA Kits include pre-coated strip plates, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents, and a detailed procedure.
Range	0.312-10 ng/ml
Sensitivity	0.095 ng/ml
Tested Applications	ELISA
Alternative Names	c-kit
Brand	Diaclone
Storage	Store at +4°C
Cross Reactivity	No cross reactivity with other human soluble receptors
Specificity	Recognizes both natural and recombinant human soluble c-kit
Sample Types	Serum, Plasma, Cell culture supernatant