

Product Datasheet

Parainfluenza Virus type 1 (PIV-1) Antibody HM226

Product Name	Parainfluenza Virus type 1 (PIV-1) Antibody
Catalog Number	HM226
Species/Host	Mouse
Tested Applications	ELISA, IFA
Brand	EastCoast Bio
Form/Appearance	Protein A Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG2a
Clonality	Monoclonal
Purity	>90%
Buffer	10 mM Phosphate Buffered Saline, pH 7.2
Specificity	Parainfluenza Virus Type 1 & 3 fusion protein
Matched Pairs	HM227
Preservatives	0.1% Sodium Azide