

# Product Datasheet

## Parainfluenza Virus type 2 (PIV-2) Antibody HM228

|                     |                                                                           |
|---------------------|---------------------------------------------------------------------------|
| Product Name        | Parainfluenza Virus type 2 (PIV-2) Antibody                               |
| Catalog Number      | HM228                                                                     |
| Species/Host        | Mouse                                                                     |
| Tested Applications | ELISA, IFA                                                                |
| Brand               | EastCoast Bio                                                             |
| Form/Appearance     | Protein A Purified Monoclonal Antibody                                    |
| Storage             | Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing. |
| Isotype             | IgG2a                                                                     |
| Clonality           | Monoclonal                                                                |
| Purity              | >90%                                                                      |
| Buffer              | 10 mM Phosphate Buffered Saline, pH 7.2                                   |
| Specificity         | Parainfluenza Virus Type 2                                                |
| Matched Pairs       | HM228 (self-pairs)                                                        |
| Preservatives       | 0.1% Sodium Azide                                                         |