

## Product Datasheet

### Troponin I (24-40) Antibody HM233

<b>Product Name</b>	Troponin I (24-40) Antibody
<b>Catalog Number</b>	HM233
<b>Description</b>	Protein A Purified and Filtered Monoclonal Antibody
<b>Species/Host</b>	Mouse
<b>Tested Applications</b>	ELISA, LF, TIA
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Protein A Purified and Filtered Monoclonal Antibody
<b>Storage</b>	+2-8 °C
<b>Isotype</b>	IgG1, kappa
<b>Clonality</b>	Monoclonal
<b>Purity</b>	≥ 95 %
<b>Buffer</b>	20 mM Sodium Phosphate, pH 7.2 with 150 mM NaCl
<b>Specificity</b>	Troponin I (epitope 24-40 a.a.)
<b>Matched Pairs</b>	HM456 (detection)