

## Product Datasheet

### Benzodiazepine (Oxazepam) Antibody HM237

<b>Product Name</b>	Benzodiazepine (Oxazepam) Antibody
<b>Catalog Number</b>	HM237
<b>Description</b>	Protein A Purified Monoclonal Antibody
<b>Tested Applications</b>	ELISA, LF
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Protein A purified monoclonal antibody
<b>Storage</b>	Short term: 2-8 °C. Long term: -20 °C. Avoid repeated freezing and thawing.
<b>Note</b>	Meets 300 ng/ml cut off for lateral flow
<b>Isotype</b>	Mouse IgG2b
<b>Clonality</b>	Monoclonal
<b>Purity</b>	≥ 95 %
<b>Buffer</b>	20 mM phosphate-buffered saline, pH 7.4 ± 0.2, with 50 mM NaCl, 0.1% sodium azide as a preservative.
<b>Specificity</b>	Oxazepam, Diazepam, and other Benzodiazepines
<b>Matched Pairs</b>	LA260 (capture)