

Product Datasheet

Benzodiazepine (Oxazepam) Antibody HM237

Product Name	Benzodiazepine (Oxazepam) Antibody
Catalog Number	HM237
Description	Protein A Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	ELISA, LF
Brand	EastCoast Bio
Form/Appearance	Protein A Purified Monoclonal Antibody
Storage	Short Term (≤ 2 weeks): +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Note	Meets 300 ng/ml cut off for lateral flow
Isotype	IgG2b
Clonality	Monoclonal
Purity	> 95 %
Buffer	20 mM Phosphate Buffered Saline, pH 7.4 with 50 mM NaCl
Specificity	Oxazepam, Diazepam, and other Benzodiazepines
Matched Pairs	LA260 (capture)
Preservatives	0.1% Sodium Azide