

Product Datasheet

Benzodiazepine (Oxazepam) Antibody HM237

Product Name	Benzodiazepine (Oxazepam) Antibody
Catalog Number	HM237
Species/Host	Mouse
Tested Applications	LF, ELISA
Brand	EastCoast Bio
Form/Appearance	Protein A Purified Monoclonal Antibody
Storage	Short Term (≤ 2 weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Note	Meets 300 ng/ml cut off for lateral flow
Isotype	IgG2b
Clonality	Monoclonal
Purity	>95%
Buffer	10 mM Phosphate Buffered Saline, pH 8.0 with 50 mM NaCl
Specificity	Oxazepam, Diazepam, and other Benzodiazepines
Matched Pairs	LA260 (capture)
Preservatives	0.1% Sodium Azide