

Product Datasheet

Human Chorionic Gonadotropin β (HCG β) Antibody HM239

Product Name	Human Chorionic Gonadotropin β (HCG β) Antibody
Catalog Number	HM239
Description	Ion Exchange Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	ELISA, LF
Brand	EastCoast Bio
Form/Appearance	Ion Exchange Purified Monoclonal Antibody
Storage	Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Note	Meets 10 mIU/ml cut off
Isotype	IgG1
Clonality	Monoclonal
Purity	> 95 %
Buffer	20 mM Phosphate Buffered Saline, pH 7.4 with 50 mM NaCl
Cross Reactivity	No reactivity with HCG α subunit or LH, TSH, FSH α and β subunits
Specificity	HCG β subunit
Matched Pairs	HM240 (capture)
Preservatives	0.1% Sodium Azide