

Product Datasheet

Luteinizing Hormone α (LH- α) Antibody HM242

Product Name	Luteinizing Hormone α (LH- α) Antibody
Catalog Number	HM242
Species/Host	Mouse
Tested Applications	LF, ELISA
Brand	EastCoast Bio
Form/Appearance	Ion Exchange Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Note	Meets 30 mIU/ml cut off
Isotype	IgG1
Clonality	Monoclonal
Purity	>95%
Buffer	10 mM Phosphate Buffered Saline, pH 8.0 with 50 mM NaCl
Cross Reactivity	No reactivity with LH, HCG, TSH, FSH β subunits
Specificity	LH α
Matched Pairs	HM241 (detection)
Preservatives	0.1% Sodium Azide