

Product Datasheet

Luteinizing Hormone α (LH- α) Antibody HM242

Product Name	Luteinizing Hormone α (LH- α) Antibody
Catalog Number	HM242
Description	Ion Exchange Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	ELISA, LF
Alternative Names	LH, Lutropin, Interstitial Cell-Stimulating Hormone
Brand	EastCoast Bio
Form/Appearance	Ion Exchange Purified Monoclonal Antibody
Storage	Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Purity	> 95 %
Buffer	10 mM Phosphate Buffered Saline, pH 8.0 with 50 mM NaCl
Cross Reactivity	No reactivity with LH, HCG, TSH, FSH β subunits
Specificity	LH α
Matched Pairs	HM241 (detection)
Preservatives	0.1% Sodium Azide