

Product Datasheet

Luteinizing Hormone α (LH- α) Antibody HM242

Product Name	Luteinizing Hormone α (LH- α) Antibody
Catalog Number	HM242
Description	Ion Exchange Purified Monoclonal Antibody
Tested Applications	ELISA, LF
Brand	EastCoast Bio
Form/Appearance	Ion exchange purified monoclonal antibody
Storage	Short term: 2-8 °C. Long term: -20 °C. Avoid repeated freezing and thawing.
Isotype	Mouse IgG1
Clonality	Monoclonal
Purity	≥ 95 %
Buffer	20 mM phosphate-buffered saline, pH 7.4 \pm 0.2, with 50 mM NaCl, 0.1% sodium azide as a preservative.
Cross Reactivity	No reactivity with LH, HCG, TSH, FSH β subunits
Specificity	LH α
Matched Pairs	HM241 (detection)