

Product Datasheet

Procalcitonin (PCT) Antibody HM243

Product Name	Procalcitonin (PCT) Antibody
Catalog Number	HM243
Description	Protein A Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	ELISA, LF
Immunogen	Recombinant procalcitonin
Brand	EastCoast Bio
Form/Appearance	Protein A Purified Monoclonal Antibody
Storage	Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Purity	> 90 %
Buffer	10 mM Phosphate Buffered Saline, pH 8.0 with 50 mM NaCl
Specificity	Human procalcitonin
Matched Pairs	HM244 (capture)
Preservatives	0.1% Sodium Azide